

MULTIPLE STIRRERS

Designed for reliable mass screenings as well as chemical, biotechnical and medical applications of all kinds, developed for the hard daily lab use and uncompromising continuous operation.

2mag MIX 15

Ultra-flat magnetic stirrer with 15 stirring positions, up to 3,000 ml stirring volume, suitable for 15 x 250 ml beaker glasses.

Powerful 15 stirring positions, 100% maintenance-free and wear-free by inductive 2mag-Magnetic-Drive concept, extremely wide speed range of 100 to 2,000 rpm, 100% synchronized speed, jerk-free stirring even at low speeds, 4-step power settings of the stirrer, high power setting for viscous liquids and reduced power setting for long-term operation without any heating effects caused by the stirrer.

Clear digital display for settings of stirrer speed and power, SoftStart procedure for reliable catching/centering and safe acceleration of the stirring bar, QuickSet for quick setting of start- and maximum speed, automatic memory function to save the last settings, outstanding long-life cycle of the hermetically sealed drive system, resistant stainless steel housing for easy and fast cleaning.

3 years warranty on material and manufacturing, developed and made in Germany.



magnetic

- 🌀 New developed, maintenance- and wear-free magnetic stirrer
- 🌀 For different vessel sizes
- 🌀 Outstanding long-life quality
- 🌀 3 years warranty

2mag MIX 6

Identical to MIX 15, but with 6 stirring positions and suitable for 6 x 1,000 ml Erlenmeyer flasks, stirring positions distance 130 mm.

magnetic^emotion



e

motion

Customized versions and modified unit sizes on request.

Technical data

Type	MIX 15	MIX 6
Stirring points	15	6
Stirring point distance	65 mm	130 mm
Stirring volume / point	1 - 3,000 ml	
Speed range	100 - 2,000 rpm	
Stirring power (max.)	20 W	
Power setting	25/50/75/100% (4-steps)	
Material housing	Stainless steel	
Operation conditions	-10 up to +50 °C (at 80% humidity)	
Electrical data	100-240 V / 50-60 Hz / 1.5 A	
Dimensions (WxDxH)	245 x 425 x 38 mm	
Protection category	IP 64	
Weight (gross)	approx. 8.8 kg	
Order no.	30015	30006

3 years warranty on material and manufacturing. We reserve the right to make technical alterations. We do not assume liability for errors in printing or typing.