
MP Biomedicals

Mycoplasma Removal Agent


One Call

One Source



Mycoplasma infected cell culture

100% removal within only one week



MRA treated cell culture

A World of Cell Culture

Eliminate Mycoplasma contamination and keep your Cells safe!

Broad spectrum activity

Ready-to-use

Reliable

Efficient

www.mpbio.com

MP Biomedicals Europe, Tel: 00800 7777 9999 • email: info.europe@mpbio.com



Mycoplasma Removal Agent (MRA)

The best solution to rescue your cells from mycoplasma contamination!

Working with human cell lines requires perfectly aseptic conditions, and often cultures get contaminated by different pathogens. Presence of mycoplasma can literally alter all cellular processes, hence the need for 100% removal. MRA is a quinolone-based antibiotic that specifically inhibits Mycoplasma DNA gyrase, an essential enzyme for microorganism replication.

- Quick: 7 days are enough to eliminate a wide variety of mycoplasma strains.
- **Efficient:** 100% removal after treatment.
- **Specific:** Inhibits mycoplasma DNA gyrase, MRA's cytotoxicity is null.
- Effortless: Simply add MRA to your standard cell culture media.

Cat. No.	Description	Pack Size
093050044	MRA – Mycoplasma Removal Agent	5 mL

For many years, MP Biomedicals has been the leading supplier of MRA. This reagent has been trusted and widely used by researchers worldwide for its excellent performance.

References:

Strand, D. et al., 2014, Am. J. Clin. Exp. Urol. 2, 239–248
Sun, G. et al., 2014, Int. J. Oncol. 45, 1091–8
Guth, A.M. et al., 2014, Vet. Immunol. Immunopathol. 161, 132–40
Pincus, S.H. et al., 2012, Vaccine 30, 4849–58

Bürdek, M. et al., 2010, J. Transl. Med. 8, 90
Habib, A. et al., 2013, Am. J. Transl. Res. 5, 634–642
Boyko-Fabian, M. et al., 2014, PLoS One 9, e98867
Sun, G.-G. et al., 2014, World J. Gastroenterol. 20, 4001–10

Also available at MP Biomedicals:

MP Biomedicals is a global leader in providing innovative tools across the **Isolate**, **Culture** and **Detect** (Immuno-solutions) workflow over 50 years.

Our goal is to provide the best quality at a competitive price combined with flexibility and customization to accelerate your research.

www.mpbio.com/immuno.php



MP Biomedicals EMEA

MP Biomedicals SARL
Parc d'Innovation, BP 50067
67402 Illkirch Cedex, France
Tel: 33 03 8867 5425
Toll Free: 800 7777 9999
Fax: 33 03 8867 1945

Email: custserv.eur@mpbio.com

MP Biomedicals LLC, HQ 3 Hutton Center Drive Suite 100, Santa Ana, California 92707 Tel: 1 949 833 2500 Toll Free: 1 800 854 0530 Fax: 1 949 859 5010 Email: sales@mpbio.com

MP Biomedicals Worldwide

MP Biomedicals APAC

MP Biomedicals
Asia Pacific Pte Ltd, HQ
2 Pioneer Place
Singapore 627885
Tel: 65 6775 0008
Fax: 65 6774 6146
Email: enquiry_ap@mpbio.com

For other offices and more information, please visit www.mpbio.com